

Ventura 4-Wheel Midsize/Standard Scooter

Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Scooter with Seating

Stadium Seat
 VENTURA418FS Standard, 18" (W) x 17.5" (D).....\$3,330.50

Captain's Seat
 VENTURA418CS Deluxe, 18" (W) x 18" (D)\$3,549.00
 VENTURA420CS Deluxe, 20" (W) x 18" (D)\$3,731.00

Batteries & Charger

12V x 20AH Batteries (pair) Included with STD
 12V x 35AH (U1) Batteries (pair)..... Included with DLX
 4A Battery Charger Included

Accessories

AB1000 Armrest Bag\$25.00
 AB1120 Backpack\$56.00

AZ2000 Scooter Cover\$72.00
 ST205-COVER 18" Seat Cover\$118.17
 ST301-COVER 20" Seat Cover\$118.17

AH1000 Crutch/Cane Holder\$53.00
 AH2000 Cylinder Tank Carrier\$70.00
 AB2100 Rear Basket\$70.00

AT1000 Scooter Trailer\$204.00

AZ0060 Cup Holder\$30.00

Specifications

HCPCS **K0807 (VENTURA418FS Only)**

Performance

Maximum Weight Capacity 350 lbs.
 Maximum Speed 5 mph
 Maximum Range 15 miles (STD) | 22 miles (DLX)
 Maximum Grade/Incline 7.5°
 Turning Radius 54"

Dimensions

Overall Length 48"
 Overall Width 25"
 Floor to Seat Height 19.5"-22.5"
 Ground Clearance 2.5"

Power

Motor 24V x 350W x 5100rpm
 Controller PG S-Drive 70A
 Brakes Electromagnetic
 Freewheel Mode Yes

Tires (Flat-free)

Front 10" x 3"
 Rear 10" x 3"

Weights

Base Weight 99 lbs.
 Battery Weight 28 lbs. (STD) | 46 lbs. (DLX)
 Seat Weight 26 lbs. (STD) | 40 lbs. (DLX)
 Total Weight 153 lbs. (STD) | 185 lbs. (DLX)

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.