

Date: _____ Account Number: _____

Company Name: _____ Contact: _____

Address: _____

_____ Phone Number: _____

PRECONFIGURED AND READY TO SHIP VIA DRIVE MEDICAL RAPID DELIVERY PROGRAM - RDP

HCPSCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims.
Referenced HCPSCS codes are also not a guarantee of payment.

WALLABY WHEELCHAIR

| ITEM # | DESCRIPTION | HCPSCS CODE | MSRP |
|--------------------------------------|---|-------------|----------|
| <input type="checkbox"/> WB1200-2GFR | Wallaby, 12" Pediatric Folding Wheelchair, Fire Truck Red | E1238 | \$550.00 |
| <input type="checkbox"/> WB1400-2GJB | Wallaby, 14" Pediatric Folding Wheelchair, Jet Fighter Blue | E1238 | \$550.00 |

Includes 2 point seat belt, flip back desk arms, swingaway footrests, and calf strap.

WALLABY WHEELCHAIR ACCESSORIES

| | | | |
|----------------------------------|--------------------------|-------|----------|
| <input type="checkbox"/> WB 8021 | Headrest Extension | E0966 | \$96.00 |
| <input type="checkbox"/> WB ELR | Elevating Legrests, Pair | E0990 | \$160.00 |
| <input type="checkbox"/> STDS819 | Anti Tippers, Pair | E0971 | \$80.00 |

*The Wallaby is crash tested and meets the ANSI/RESNA WC19 crash test standards when used together with the headrest extension for user's between 50 and 80 pounds.

WALLABY WHEELCHAIR SPECIFICATIONS

| | WB1200-2GFR | WB1400-2GJB |
|--|-------------------|-------------------|
| Overall Width (open closed) | 20" 12" | 22" 12" |
| Seat Dimensions | 12" (W) x 12" (D) | 14" (W) x 14" (D) |
| Seat-to-Floor (front) | 19" | 19" |
| Seat-to-Footrest | 13" - 16" | 13" - 16" |
| Back Upholstery Height | 14" | 14" |
| Push Handle Height | 33" - 39" | 33" - 39" |
| Armrest Height | 7" | 7" |
| Overall Length with Footrests | 37" | 39" |
| Overall Length without Footrests | 28.5" | 30.5" |
| Product Weight with Footrests | 37 lbs. | 37 lbs. |
| Product Weight without Footrests | 33 lbs. | 33 lbs. |
| Weight Capacity | 150 lbs. | 150 lbs. |
| Weight Capacity when Used on Bus Transport | 50 - 80 lbs. | 50 - 80 lbs. |