


770-537 AQUASENSE ADJUSTABLE BATH BENCH WITH BACK, white



Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	A11-452	Plastic Plug for dia:22.2mm, Aqua-Blue - 4ea / set	1	1	ST
2.	B11-832	Hardware Set, Aqua-Blue (5 Triangular Knobs M6x34mm & 4 Triangular Knobs M6x57mm)	1	1	ST
3.	B11-811	Push Button dia:8mm & Spring - 4ea / set	1	1	ST
4.	A11-212	Anti-Rattle for dia:28.6mm, Aqua-Blue - 4ea / set	1	1	ST
5.	B11-012	Extension Leg, Aqua-Blue (Aluminum Tube dia:28.6x250mm, Anti-Rattle & Rubber Tip) - 4ea / set	1	1	ST
6.	A11-482	Rubber Tip for dia:28.6mm, Aqua-Blue - 4ea / set	1	1	ST
7.	A11-426	Backrest, Blow-Molded Plastic, White (410x365x25mm)	1	1	EA
8.	A11-406	Seat, Blow-Molded Plastic, White (500x305x40mm)	1	1	EA
Not Shown	B11-852	Hardware Set for Bath Bench with Back (5 Cross-Cup Screws M6x34mm, 4 Cross-Cup Screws M6x59mm, 9 Curved Washers) - Accessory	1	1	ST