


770-120 FOLDING WALKER WITH 5" WHEELS ADULT, SILVER



Item #	AMG Part #	Item Description	Qty Per Unit	Sold by	UOM
1.	A03-402YE	Push-Button Cover, Yellow	1	1	EA
2.	770-168	Handle, Grey	2	1	EA
3.	B03-841	Hardware for Folding Bar (Bolt M6x66mm, 3 Grey Curved Washers, 1 Flat Washer & 1 Hex Nut)	2	1	ST
4.	A03-210	Anti-Rattle, Grey	8	1	EA
5.	B03-811	Push-Button & Spring	4	1	ST
6.	780-118	Walker Glides, Grey	1	1	PR
7.	770-162	5" Wheel X-Long Extension for 1" Tube, Silver	1	1	PR
8.	B03-851GY	End Plug Set, Grey (3/4" Plug & 7/8" Plug)	1	1	ST