



# Premium Raised Toilet Seat with Removable Arms, Standard Seat

Item # 12402



STANDARD ROUND




<b>A</b>	Overall Seat Width	13.75"
<b>B</b>	Inside Seat Width	9.75"
<b>C</b>	Overall Seat Depth/Length	16.5"
<b>D</b>	Inside Seat Depth/Length	11.25"

<b>E</b>	Overall Seat Height	3.5"
<b>F</b>	Overall Width (with Handles)	21.75"
<b>G</b>	Overall Height (with Handles)	9.25"
<b>H</b>	Mounting Hole Width	5.5"

# Premium Raised Toilet Seat with Removable Arms, Elongated Seat

Item # 12403




<b>A</b>	Overall Seat Width	14"
<b>B</b>	Inside Seat Width	10"
<b>C</b>	Overall Seat Depth/Length	18.5"
<b>D</b>	Inside Seat Depth/Length	13"

<b>E</b>	Overall Seat Height	3.5"
<b>F</b>	Overall Width (with Handles)	21"
<b>G</b>	Overall Height (with Handles)	10"
<b>H</b>	Mounting Hole Width	5.5"

# Premium Plastic, Raised, Regular/Elongated Toilet Seat with Lock

Item # 12013




<b>A</b>	Overall Seat Width	15.5"
<b>B</b>	Inside Seat Width	8"
<b>C</b>	Overall Seat Depth/Length	16.5"
<b>D</b>	Inside Seat Depth/Length	9"

<b>E</b>	Overall Seat Height	5"
<b>F</b>	Overall Width (with Handles)	22.5"
<b>G</b>	Overall Height (with Handles)	9"
<b>H</b>	Seat Lock Range	11.5"-15.5"

# Premium Plastic Raised, Regular/Elongated Toilet Seat, with Lock

Item # 12014



<b>A</b>	Overall Seat Width	17"
<b>B</b>	Inside Seat Width	8"
<b>C</b>	Overall Seat Depth/Length	16.5"

<b>D</b>	Inside Seat Depth/Length	9"
<b>E</b>	Overall Seat Height	4.5"
<b>F</b>	Seat Lock Range	11.5-15.5"