



Raised Toilet Seat with Lock, 2"

Item # 12062



STANDARD ROUND




A	Overall Seat Width	14"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	15.5"


D	Inside Seat Depth/Length	10"
E	Overall Seat Height	2.5"
F	Mounting Hole Width	N/A

Raised Toilet Seat with Lock, 6"

Item # 12066



STANDARD ROUND




A	Overall Seat Width	14"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	15.5"

D	Inside Seat Depth/Length	10"
E	Overall Seat Height	6"
F	Mounting Hole Width	N/A

Raised Toilet Seat with Lock, 4"

Item # RTL12064




A	Overall Seat Width	14"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	15.5"

D	Inside Seat Depth/Length	10"
E	Overall Seat Height	4"
F	Mounting Hole Width	N/A

Raised Toilet Seat with Lock and Lid, 2"

Item # 12063




A	Overall Seat Width	14"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	16"

D	Inside Seat Depth/Length	10.5"
E	Overall Seat Height	2.5"
F	Mounting Hole Width	N/A

Raised Toilet Seat with Lock and Lid, 4"

Item # 12065




A	Overall Seat Width	14"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	16"

D	Inside Seat Depth/Length	10"
E	Overall Seat Height	4"
F	Mounting Hole Width	N/A

Raised Toilet Seat with Lock and Lid, 6"

Item # 12067



A	Overall Seat Width	14"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	16"

D	Inside Seat Depth/Length	10"
E	Overall Seat Height	6"
F	Mounting Hole Width	N/A