


Hinged Toilet Seat Riser, Standard

Item # RTL12607




A	Overall Seat Width	13.5"
B	Inside Seat Width	9"
C	Overall Seat Depth/Length	16.5"

D	Inside Seat Depth/Length	10.75"
E	Overall Seat Height	3.5"
F	Mounting Hole Width	5.5"

Hinged Toilet Seat Riser, Elongated

Item # RTL12608



A	Overall Seat Width	14"
B	Inside Seat Width	10"
C	Overall Seat Depth/Length	18.5"

D	Inside Seat Depth/Length	13.5"
E	Overall Seat Height	3.5"
F	Mounting Hole Width	5.5"