

Drive Titan LTE Portable Powerchair Order Form/Price List



	Address:	
	Phone #:	Fax #:
	Email:	

Account:		Phone #:
Contact Person:	P.O. #:	Date:
Bill To:	Ship To:	

PORTABLE POWERCHAIR		
Folding Seat		
<input type="checkbox"/>	TITANLTE-18FS	18" (W) X 18" (D)
BATTERIES AND CHARGER		
	12v x 20AH (pair)	Included
	3A Battery Charger (offboard)	Included
ACCESSORIES		
<input type="checkbox"/>	SF8010 Oxygen Tank Carrier	
<input type="checkbox"/>	SF8020 Rear Basket	
<input type="checkbox"/>	SF8030 Cane/Crutch Holder	
<input type="checkbox"/>	AE2500 ELR Bracket, Hemi Spacing	
<input type="checkbox"/>	LK3JELR ELR Set	

SPECIFICATIONS	
HCPCS	
TITANLTE-18FS	K0821
PERFORMANCE	
Maximum Weight Capacity	300 lbs.
Maximum Speed	4.2 mph
Maximum Range	8 miles
Maximum Grade/Incline	6°
Turning Radius	27"
DIMENSIONS	
Seat Width	18"
Seat Depth	18"
Floor to Seat Height	20 - 22"
Ground Clearance	2.5
POWER	
Motor	200W x 3900 rpm
Controller	Dynamic Linx 40A
Brakes	Electromagnetic
Freewheel Mode	Yes
TIRES (Flat-free)	
Drive Wheels	9" x 3"
WEIGHTS	
Base Weight w/Out Battery	70 lbs.
Battery Weight	22 lbs.
Seat Weight	30 lbs.
Heaviest Piece Weight	39 lbs.
Total Weight	122 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.