

## Disposable Container with Internal Filter



**Assembled Collection Container Kit 7305D-633** - Includes internal filter cartridge, 800 ml container and lid, splash guard, and 4.5" connection and 6' patient tubing (1 ea)

**Assembled Collection Container Kit 7305D-633-25** - Includes internal filter cartridge, 800 ml container and lid, splash guard, and 4.5" connection and 6' patient tubing (25 pk)



**Unassembled Collection Container 7305D-632** - Includes internal filter cartridge, 800 ml container and lid, splash guard, and 4.5" connection tubing (48 pk)



**Elbow with Integrated Filter 7305D-635** (12/pk)



**6' Patient Tubing SUCP TUBING 72** (1/ea)

## Disposable Container with External Filter



**Assembled Collection Container Kit 22330** - Includes 800 ml disposable container & lid, external bacteria filter, elbow, 10" connection and 6' patient tubing package (1/ea)

**Assembled Collection Container Kit 22330-12** - 800 ml disposable container & lid, external bacteria filter, elbow, 10" connection and 6' patient tubing package (12/pk)



**Disposable Container and Lid only 610-12B** (12/pk)  
610-48BP (48/pk)



**6' Patient Tubing SUCP TUBING 72** (1/ea)



**Universal Suction Tubing and Filter Kit 18600-KITN** - (filter, elbow, 10" connection and 6' patient tubing package) (1/ea)



**3.375" Connection Tubing 7305D-612** (6/pk)



**Filter 18600-FILTER** (3/pk)



**Filter 7305D-608** (12/pk)



**Elbow 7305D-619** (6/pk)



**Elbow, Filter, 3.375" Connection Tubing Kit 7305D-614** (1/ea)

**NOTE** - The 18600 Drive Heavy-Duty Suction Unit requires the 10" tubing connection filter included in the 22330 container series only. Any container series may be used with the 7305 and 7314 DeVilbiss Vacu-Aide® Suction Units.

## Accessories and Replacement Parts



**7305D Carrying Case 7305D-606**



**12V DC Power Cord 7304D-619**



**AC to DC Power Adapter/Charger 7314P-613**



**Power Cord for USA DV51D-606**



**Vacuum Gauge Test Fixture 7304D-636**



**Connection Tubing SUCP TUBING 10**

**Battery with Foam for 7314P-D 7314P-614**

**Battery for 7305P-D 7305P-614**

**Battery Retrofit Kit for 7305D-D 7305P-620**



**7314D Carrying Case 7314D-606**