



**PreserveTech™**  
Stays Cleaner • Lasts Longer

# Universal Raised Toilet Seat

Introducing the **PreserveTech™** Universal Raised Toilet Seat, our first-ever truly universal raised toilet seat accommodates three different seat heights (2, 4 or 6") by maneuvering just two components. This seat can fit on most standard and elongated toilets and features two hygienic cutaways for easy cleaning. The pivoting armrests can be raised or removed, depending on user needs and bathroom space. This item also includes a lid with a recessed grip for effortless opening.

Drive **PreserveTech™** products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.

\*Antimicrobial properties are built-in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.



HOW TO CHOOSE A RAISED TOILET SEAT



*Scan me*

ITEM #	PRODUCT DESCRIPTION	UPC	UOM
RTL12C002-WH	PreserveTech™ Universal Raised Toilet Seat	822383006147/ea / 50822383006142/cv	1/ea, 1/cv

### PRODUCT SPECS

Product Dimensions	22.00" (W) x 18.50" (L) x 15.00" (H)
Product Weight	9.30 lb.
Color	White
Weight Capacity	350 lb.
Warranty	Limited Lifetime**

### BOX SPECS

Retail Box Dimensions	8.66" (W) x 20.28" (L) x 18.50" (H)
Retail Box Weight	11.07 lb.

\*\*Refer to instruction guide or website for warranty details.