


## TITAN (Slide out battery tray after serial #(21)2W1808202659)



- |  |   |
|--|---|
| 1. Footplate & Attaching Bolt- TITAN-01<br>Footplate Depth Adj. Bolt- TITAN-14                   | 25. Captain's Seat Assembly-<br>16X16X18- TITAN-29<br>16X18X18- TITAN-30<br>18X16X18- TITAN-31<br>18X18X16.5- TITAN-19<br>18X18X18- TITAN-20<br>20X18X18- TITAN-21<br>20X20X18- TITAN-32<br>22X20X19.5- TITAN-33<br>Pan Seat-<br>18"- ST501-181818<br>20"- ST502-202020 |
| 2. Frame- C01-097-00100  | 26. Headrest- ST-HEADREST   |
| 3. Anti-Tipper Wheel- TITAN-41   | 27. Seat Belt- STDS855-PW   |
| 4. Battery Tray Bracket & Hdw.- TITAN-43   | 28. Armrest Assembly Left- LDR410011L<br>Armrest Assembly Right- LDR410011R   |
| 5. Battery Strap (ea)- C31-005-00400   | 29. Arm Pad- LDR410110  |
| 6. Battery (ea)(2 needed)- LRA402208   | 30. Armrest Height & Width Knob- S31052   |
| 7. Battery Wire Harness- TITAN-12<br>Circuit Breaker- LRA402307                                  | 31. Seat Plate- C07-072-01001<br>Seat Plate Accessory Knob- S31052  |
| 8. Motor & Gear Box Left- TITAN-10<br>Motor & Gear Box Right- TITAN-09<br>Axle Key Way- TITAN-40 | 32. Seat Post- TITAN-39 [1]<br>Seat Post Bolt- TITAN-06   |
| 9. Crutch Handle Assembly- TITAN-26<br>Crutch Handle Tip- TITAN-27                               | 33. Joystick Bracket- TITAN-24<br>Joystick Bracket Knob- S31052   |
| 10. Drive Wheel Complete- TITAN-11   | 34. Joystick Knob- LDM400001-04<br>Joystick Mounting Screws (pr)- LDM400001-01  |
| 11. Drive Wheel Nut & Washer- TITAN-28   | 35. Joystick- LINX-JOYSTICK<br>Joystick/Controller Cable- GLM-BUS150-A<br>Joystick Tip- DC GPL55470   |
| 12. Hub Cap- LRM412202-HC  | 36. Charger- LRA402209  |
| 13. Caster- TITAN-03<br>Caster Bearing- STDS1011   |   |
| 14. Caster Axle- TITAN-04  |   |
| 15. Fork- LDM400013  |   |
| 16. Dust Cap- LDR4110  |   |
| 17. Fork Bearing Set- LRB402348  |   |
| 18. Main Wire Harness- C18-097-00100   |   |
| 19. Controller- DLX-PM40-A-TITAN   |   |
| 20. Front Shroud- C08-097-00100  |   |
| 21. Side Shroud Left- C08-097-00700<br>Side Shroud Right- C08-097-00600                          |   |
| 22. Rear Shroud Cover- C08-097-00200   |   |
| 23. Rear Bumper- C08-097-00500   |   |
| 24. Decoration Set Blue- TITAN-18N<br>Decoration Set Red- TITAN-17N                              |   |

[1] No longer available in chrome once stock is gone.