

Scout 4-Wheel Transport Scooter Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Scooter with Seating

Folding Seat

- SFSCOUT4.....16.5" (W) x 16" (D)\$1,820.00
HCPCS: K0800, Includes interchangeable color shrouds in Blue & Red.
- SFSCOUT4-EXT.....16.5" (W) x 16" (D)\$1,929.20
HCPCS: K0800, Includes interchangeable color shrouds in Blue & Red.

Batteries & Charger

- 12V x 12AH (pair)Included
- 12V x 21AH (pair)Included (EXT Model)
- 2A Battery ChargerIncluded

Accessories

- SF8010 Cylinder Tank Carrier (HCPCS E2208)\$150.91
- SF8020 Rear Basket\$191.89
- SF8030 Crutch/Cane Holder (HCPCS E2207)\$101.23
- SF8050 Sun Shade\$661.20

Specifications

HCPCS

SFSCOUT4	K0800
SFSCOUT4-EXT	K0800

Performance

Maximum Weight Capacity	300 lbs.
Maximum Speed	4.25 mph
Maximum Range	9 miles
Maximum Range (EXT Model)	15 miles
Maximum Grade/Incline	6°
Turning Radius	53.75"

Dimensions

Overall Length	42.5"
Overall Width	21.25"
Floor to Seat Height	17 - 19"
Ground Clearance	2.5"

Power

Motor	24V x 270W x 4700rpm
Controller	45A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires (Flat-Free)

Front	8" x 2"
Rear	8" x 2"

Weights

Base Weight	58 lbs.
Battery Weight	20 lbs.
Seat Weight	16 lbs.
Heaviest Piece Weight	39 lbs. (Front Section)
Total Weight	94 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.