



Premium Raised Toilet Seat with Removable Arms, Standard Seat

Item # 12402



STANDARD ROUND




A	Overall Seat Width	13.75"
B	Inside Seat Width	9.75"
C	Overall Seat Depth/Length	16.5"
D	Inside Seat Depth/Length	11.25"

E	Overall Seat Height	3.5"
F	Overall Width (with Handles)	21.75"
G	Overall Height (with Handles)	9.25"
H	Mounting Hole Width	5.5"

Premium Raised Toilet Seat with Removable Arms, Elongated Seat

Item # 12403




A	Overall Seat Width	14"
B	Inside Seat Width	10"
C	Overall Seat Depth/Length	18.5"
D	Inside Seat Depth/Length	13"

E	Overall Seat Height	3.5"
F	Overall Width (with Handles)	21"
G	Overall Height (with Handles)	10"
H	Mounting Hole Width	5.5"

Premium Plastic, Raised, Regular/Elongated Toilet Seat with Lock

Item # 12013




A	Overall Seat Width	15.5"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	16.5"
D	Inside Seat Depth/Length	9"

E	Overall Seat Height	5"
F	Overall Width (with Handles)	22.5"
G	Overall Height (with Handles)	9"
H	Seat Lock Range	11.5"-15.5"

Premium Plastic Raised, Regular/Elongated Toilet Seat, with Lock

Item # 12014



A	Overall Seat Width	17"
B	Inside Seat Width	8"
C	Overall Seat Depth/Length	16.5"

D	Inside Seat Depth/Length	9"
E	Overall Seat Height	4.5"
F	Seat Lock Range	11.5-15.5"