

# Phoenix HD 4-Wheel Transport Scooter Order Form/Price List



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>	<b>P.O. #</b>	<b>Date:</b>	
<b>Bill to:</b>		<b>Ship to:</b>	

### Scooter with Seating

#### Folding Seat

- PHOENIXHD4 .....17.5" (W) x 17" (D) .....\$2,730.00  
*HCPCS: K0801, Includes interchangeable color shrouds in Blue & Red.*
- PHOENIXHD4-20.....20" (W) x 17" (D) .....\$2,857.40  
*HCPCS: K0801, Includes interchangeable color shrouds in Blue & Red.*

### Batteries & Charger

- 12V x 20AH (pair) ..... Included
- 1.5A Battery Charger ..... Included

### Accessories

- SF8010 Cylinder Tank Carrier (HCPCS E2208) .....\$150.91
- SF8020 Rear Basket .....\$191.89
- SF8030 Crutch/Cane Holder (HCPCS E2207) .....\$101.23
- SF8050 Sun Shade .....\$661.20

### Specifications

#### HCPCS

PHOENIXHD4	K0801
PHOENIXHD4-20	K0801

#### Performance

Maximum Weight Capacity	350 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	54"

#### Dimensions

Overall Length	41.5"
Overall Width	22.5"
Floor to Seat Height	20 - 22"
Ground Clearance	3.75"

#### Power

Motor	24V x 350W x 3600rpm
Controller	70A
Brakes	Electromagnetic
Freewheel Mode	Yes

#### Tires (Flat-Free)

Front	9" x 3"
Rear	9" x 3"

#### Weights

Base Weight	77 lbs.
Battery Weight	31 lbs.
Seat Weight	21 lbs.
Heaviest Piece Weight	46 lbs. (Front Section)
Total Weight	129 lbs.

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.*