


## TITAN AXS



- |  |  |
|--|--|
| 1. Shroud- C08-084-00100   | 27. Footrest Mtg Bracket Assembly- TITANMID-04 |
| 2. Shroud Knob- VENT-86  | 28. Main Frame- TITANMID-05                    |
| 3. Decoration Set- S08-084-00118   | 29. Frame Bushing Assembly- TITANMID-06        |
| 4. Front Cover & Hardware- TITANMID-01   | 30. Side Frame Assembly Left- TITANMID-07      |
| 5. Motor, Gear Brake Assy. Left- C09-084-00201                                   | 31. Side Frame Assembly Right- TITANMID-08     |
| 6. Motor, Gear Brake Assy. Right- C09-084-00101                                  | 32. Seat Post *- TITAN-39                      |
| 7. Joystick w/o Extension Cable- LINX-JOYSTICK<br>Joystick Tip Linx- DC GPL55470 | 33. Seat Post Bolt Assembly- TITAN-06          |
| 8. Joystick Arm- LDM400001-04<br>Joystick Mtg. Screw- LDM400001-01               | 34. Seat Plate- C07-072-01001*                 |
| 9. Joystick Mtg. Bracket- TITAN-24   | 35. Captain Seat 16 x 16- TITANAXS-01*         |
| 10. Joystick Knob- LDM400001-07  | Captain Seat 16 x 18- TITANAXS-02*             |
| 11. Joystick Extension Cable- GLM-EXT100-A                                       | Captain Seat 18 x 16- TITANAXS-03*             |
| 12. Joystick Velcro Cable Tie- C34-072-00100                                     | Captain Seat 18 x 18- TITANAXS-04*             |
| 13. Controller- DLX-PM40-A-AXS<br>Controller to Shroud Cable- GLM-EXT064-A       | Captain Seat 20 x 18- TITANAXS-05*             |
| 14. Wire Harness- TITAN-12   | Captain Seat 20 x 20- TITANAXS-06*             |
| 15. Battery 35 ah (2 needed)- LRA402208  | Captain Seat 22 x 20- TITANAXS-07*             |
| 16. Charger 3 Amp- LRA402209   | Pan Seat 18 x 18- ST501-181818-AXS             |
| 17. Caster Axle & Nut- TITAN-04  | Pan Seat 20 x 18- ST502-201820-AXS             |
| 18. Caster- TITAN-03<br>Caster Bearing- STDS1011                                 | 36. Headrest- ST-HEADREST                      |
| 19. Fork- LDM400013  | 37. Armrest Assembly Left- LDR410011L          |
| 20. Fork Bearing Set Assembly- LDB402348   | Armrest Assembly Right- LDR410011R             |
| 21. Fork Cap- LDR4110  | Arm Pad- LDR410110                             |
| 22. Drive Wheel- TITAN-11  | Arm Pad Screw (pr)- LDR420018                  |
| 23. Axle Nut- EC16018  | 38. Armrest Height Adj. Knob- S31052           |
| 24. Hub Cap- LRM412202-HC  | 39. Armrest Width Adj. Knob- S31052            |
| 25. Footrest Hanger Assembly- TITANMID-02  | 40. Seatbelt 48"- STDS850-PW                   |
| 26. Footrest Assembly- TITANMID-03   | 60"- STDS855-PW                                |
|  | Not Shown-                                     |
|  | ELR Adapter- AE2500                            |
|  | ELR - LK3JELR                                  |

**\* When current stock is gone the part will no longer be chrome plated. It will be black.**