

Scout 3-Wheel Transport Scooter Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Scooter with Seating

Folding Seat

- SFSCOUT3.....16.5" (W) x 16" (D)\$1,765.50
HCPCS: K0800, Includes interchangeable color shrouds in Blue & Red.
- SFSCOUT3-EXT.....16.5" (W) x 16" (D)\$1,874.60
HCPCS: K0800, Includes interchangeable color shrouds in Blue & Red.

Batteries & Charger

- 12V x 12AH (pair) Included
- 12V x 20AH (pair) Included (EXT Model)
- 2A Battery Charger Included

Accessories

- SF8010 Cylinder Tank Carrier (HCPCS E2208)\$150.91
- SF8020 Rear Basket\$191.89
- SF8030 Crutch/Cane Holder (HCPCS E2207)\$101.23
- SF8050 Sun Shade\$661.20

Specifications

HCPCS

SFSCOUT3	K0800
SFSCOUT3-EXT	K0800

Performance

Maximum Weight Capacity	300 lbs.
Maximum Speed	4.25 mph
Maximum Range	9 miles
Maximum Range (EXT Model)	15 miles
Maximum Grade/Incline	6°
Turning Radius	45.5"

Dimensions

Overall Length	42.5"
Overall Width	21.25"
Floor to Seat Height	17 - 19"
Ground Clearance	2.5"

Power

Motor	24V x 270W x 4700rpm
Controller	45A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires (Flat-Free)

Front	8" x 2"
Rear	8" x 2"

Weights

Base Weight	49 lbs.
Battery Weight	20 lbs.
Seat Weight	16 lbs.
Heaviest Piece Weight	29 lbs. (Front Section)
Total Weight	85 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims for reimbursement and should consult an attorney or other advisor to discuss specific situations in further detail.