

Phoenix HD 3-Wheel Transport Scooter Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Scooter with Seating

Folding Seat

- PHOENIXHD317.5" (W) x 17" (D)\$2,548.00
HCPCS: K0801, Includes interchangeable color shrouds in Blue & Red.

Batteries & Charger

- 12V x 20AH (pair) Included
- 1.5A Battery Charger Included

Accessories

- SF8010 Cylinder Tank Carrier (HCPCS E2208)\$150.91
- SF8020 Rear Basket\$191.89
- SF8030 Crutch/Cane Holder (HCPCS E2207)\$101.23
- SF8050 Sun Shade\$661.20

Specifications

HCPCS

PHOENIXHD3	K0801
PHOENIXHD3-20	K0801

Performance

Maximum Weight Capacity	350 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	43"

Dimensions

Overall Length	39"
Overall Width	22.5"
Floor to Seat Height	20 - 22"
Ground Clearance	3.75"

Power

Motor	24V x 350W x 3600rpm
Controller	70A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires (Flat-Free)

Front	9" x 3"
Rear	9" x 3"

Weights

Base Weight	63 lbs.
Battery Weight	31 lbs.
Seat Weight	21 lbs.
Heaviest Piece Weight	32 lbs. (Front Section)
Total Weight	115 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims for reimbursement and should consult an attorney or other advisor to discuss specific situations in further detail.