




PHOENIX HD3

Delta Tiller s/n (21)2W1801100007 and later.	Standard tiller before s/n (21)2W1801100007
 <ol style="list-style-type: none"> 1. Throttle- S12-034-00300 2. Speed Control- S31001 3. Charger Port- S31147 4. Key Switch- S31006 KEY- S3KEY 5. Delta Bar Cover- S35010HD-03 (Upper & Lower Cover) 6. Headlight Assembly- S35010HD-02 7. Front Tiller Cover- C08-034-01100 <ol style="list-style-type: none"> a. Battery Gauge LED- C11-076-00100 b. Headlight Switch- S31122 c. Horn Button- S31012 d. USB Assembly- S35010HD-01 e. Control Panel Assembly- C08-034-01200 8. Rear Tiller Cover Upper- C08-034-01300 Lower- C08-034-01400 	 <ol style="list-style-type: none"> 1. Throttle Old- * New- S12-069-00300 * Stock is gone use new style. 2. Speed Control- S31001 3. Charger Port- S31147 4. Key Switch- S31006 Key- S3KEY 5. Handgrip Foam- S31003 Handgrip TPR Plastic- C20-007-00200 6. Headlight Assembly- S31022 7. Front Tiller Cover- S31004GT <ol style="list-style-type: none"> a. Battery Gauge LED- C11-076-00100 LED- C11-076-00100 Analog- C11-033-00101 b. Headlight Switch- S31122 c. Horn Button- S31012 8. Lower Tiller Cover- S31209BK 9. Tiller Standard- S31002 Tiller Ext. Foot Platform- C02-034-00100



- | | |
|---|--|
| <p>10. Basket Bracket- BC35051</p> <p>11. Tiller Swivel Assy.- S350109
Tiller Adj. Knob- S31028</p> <p>12. Tiller Bearing Assy.- S31021</p> <p>13. Fork- S31036-3HD
Ext Foot Platform- C03-009-00101</p> <p>14. Fork Cover- S350127</p> <p>15. Blue Decoration- S35012B
Red Decoration- S35012R</p> <p>16. Front Wheel- S31142-3HD
Front Wheel Bearing- 9503H79021
10mm ID – 22mm OD</p> <p>17. Front Axle- S05-009-00100</p> <p>18. Tow Handle- S31011</p> <p>19. Tension Knob- S350169715</p> <p>20. T Handle Frame lock- S350169716</p> <p>21. Seat pin- S350169717</p> <p>22. Main Shroud Decoration
Blue Left- S3501BL
Blue Right- S3501BR
Red Left- S3501RL
Blue Right- S3501RR</p> <p>23. Main Shroud- S350110</p> <p>24. Controller Cover- S350113</p> <p>25. Foot Shroud- S350114
Ext Foot Platform- C08-035-01500</p> <p>26. Battery Pack- S350169707-HD
Case Only- S350169706-HD
Charger Port- S350169706-HDCP</p> <p>27. Battery (2 needed)- LRM412123</p> <p>28. Charger- LRM302111</p> <p>29. Controller (CTE)- S31008GT-HD
Controller (S-Drive)- C10-067-00200 *
* After Serial # (21)2W1801100007</p> <p>30. Rear Wheel- S31132B-HD</p> <p>31. Rear Shroud- S350111
Lifting Handle- S35011FH</p> | <p>32. Rear Shroud Decoration
Blue Left- S35011B-L
Blue Right- S35011B-R
Red Left- S35011R-L
Red Right- S35011R-R</p> <p>33. Front Frame
Non-Ext Foot Platform- NLA
Ext Foot Platform- C01-035-00800</p> <p>34. Rear Frame
Non-Ext Foot Platform- S31130-4HD
Ext Foot Platform- C01-035-00900</p> <p>35. Seat Post- S31049</p> <p>36. Wire Harness (CTE)- S31048-HD
Wire Harness (S Drive)- S18-034-00400</p> <p>37. Rear Axle Nut- S350112AN</p> <p>38. Motor/Brake/Transaxle- S31125HD
Motor/Brake- S31031-HD
Transaxle- S350112-HD
Axle Key- 29002W7011</p> <p>39. Transaxle Bracket- S350112BKT-HD</p> <p>40. Anti-tip Wheel-S31026</p> <p>41. Armrest Assy. Left- VENT-29
Right- VENT-30
Armrest Pad Left- VENT-31
Right- VENT-32</p> <p>42. Seat Plate- S31039-HD</p> <p>43. Armrest- Adj. Knob- S31052</p> <p>44. 17.5" Seat Assembly- S35045-HD
17.5" Back Cover- S35045BC-HD
17.5" Seat Cover- S35045SC-HD
20" Seat & Seat Plate- S35045HD20</p> <p>45. Frame Contact Front- S31056HD
Rear- S31057HD</p> <p>46. Basket- S35003-PL</p> <p>47. Free Wheel Knob- S31031K</p> <p>48. Mud Flap Left- S35048MFL
Right- S35048MFR</p> <p>49. Battery Contact- S31060</p> <p>50. Battery Box Charge Port- S350169706-HDCP</p> |
|---|--|