

Date: \_\_\_\_\_ Account Number: \_\_\_\_\_

Company Name: \_\_\_\_\_ Contact: \_\_\_\_\_

Address: \_\_\_\_\_

Phone Number: \_\_\_\_\_

**PRECONFIGURED AND READY TO SHIP VIA DRIVE MEDICAL RAPID DELIVERY PROGRAM - RDP**

HCPSCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims.  
Referenced HCPSCS codes are also not a guarantee of payment.

**WALLABY WHEELCHAIR**

ITEM #	DESCRIPTION	HCPSCS CODE	MSRP
<input type="checkbox"/> WB1200-2GFR	Wallaby, 12" Pediatric Folding Wheelchair, Fire Truck Red	E1238	\$550.00
<input type="checkbox"/> WB1400-2GJB	Wallaby, 14" Pediatric Folding Wheelchair, Jet Fighter Blue	E1238	\$550.00

Includes 2 point seat belt, flip back desk arms, swingaway footrests, and calf strap.

**WALLABY WHEELCHAIR ACCESSORIES**

<input type="checkbox"/> WB 8021	Headrest Extension	E0966	\$96.00
<input type="checkbox"/> WB ELR	Elevating Legrests, Pair	E0990	\$160.00
<input type="checkbox"/> STDS819	Anti Tippers, Pair	E0971	\$80.00

\*The Wallaby is crash tested and meets the ANSI/RESNA WC19 crash test standards when used together with the headrest extension for user's between 50 and 80 pounds.

**WALLABY WHEELCHAIR SPECIFICATIONS**

	WB1200-2GFR	WB1400-2GJB
Overall Width (open   closed)	20"   12"	22"   12"
Seat Dimensions	12" (W) x 12" (D)	14" (W) x 14" (D)
Seat-to-Floor (front)	19"	19"
Seat-to-Footrest	13" - 16"	13" - 16"
Back Upholstery Height	14"	14"
Push Handle Height	33" - 39"	33" - 39"
Armrest Height	7"	7"
Overall Length with Footrests	37"	39"
Overall Length without Footrests	28.5"	30.5"
Product Weight with Footrests	37 lbs.	37 lbs.
Product Weight without Footrests	33 lbs.	33 lbs.
Weight Capacity	150 lbs.	150 lbs.
Weight Capacity when Used on Bus Transport	50 - 80 lbs.	50 - 80 lbs.