

Titan Front Wheel Drive Power Wheelchair Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Powerchair with Seating

Standard Back Captain's Seat

- TITANLB18CSX1618" (W) x 18" (D)\$5,560.00
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.
- TITAN1616X1616" (W) x 16" (D)\$5,395.00
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.
- TITAN1618X1616" (W) x 18" (D)\$5,460.50
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.
- TITAN1816X1618" (W) x 16" (D)\$5,494.00
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.

High Back Captain's Seat

- TITAN18CSX1618" (W) x 18" (D)\$5,725.50
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.
- TITAN20CSX1620" (W) x 18" (D)\$5,891.50
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.
- TITAN2020X1620" (W) x 20" (D)\$5,957.00
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.
- TITAN22CSX1622" (W) x 20" (D)\$6,056.00
HCPCS: K0816, Includes interchangeable color shrouds in Blue & Red.

Standard Armrests

- ARM-L.....Included
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests
- ARM-RIncluded
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests

Controller Arm

- Adjustable Length.....Included
- AA4100 – Swing-Away Controller Arm (HCPCS E1028)\$105.00

Front Riggings

- Adjustable FootplateIncluded
- AE2500/LK3JELR – Elevating Legrest Kit\$207.48
HCPCS: E0990, Includes swing-away elevating legrest with attachment bracket

Batteries & Charger

- U1BATTERY (pair)Included
HCPCS: E2365, 12V x 35AH
- 3A Battery ChargerIncluded

Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978)Included
- AH1000 Crutch/Cane Holder (HCPCS E2207)\$120.12
- AH2000 Cylinder Tank Carrier (HCPCS E2208)\$156.52
- ST205-COVER 18" Seat Cover\$118.17
- ST301-COVER 20" Seat Cover\$118.17
- AZ0060 Cup Holder\$50.96

Specifications

HCPCS	
TITANLB18CSX16	K0816
TITAN1616X16	K0816
TITAN1618X16	K0816
TITAN1816X16	K0816
TITAN18CSX16	K0816
TITAN20CSX16	K0816
TITAN2020X16	K0816
TITAN22CSX16	K0816
Performance	
Maximum Weight Capacity	300 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	25"
Dimensions	
Overall Length	40"
Overall Width	25"
Floor to Seat Height	20.5"–23"
Ground Clearance	2.5"
Power	
Motor	Two 24V x 200W x 3900rpm
Controller	Dynamic Shark 40A
Brakes	Electromagnetic
Freewheel Mode	Yes
Tires	
Casters	6" x 2" Solid
Drive Wheels	10" x 3" Flat Free
Weights	
Base Weight	68 lbs.
Battery Weight	44 lbs.
Seat Weight	43 lbs.
Heaviest Piece Weight (Disassembled)	54 lbs.
Total Weight	155 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.