



Phoenix HD4



- | | |
|---|--|
| <ol style="list-style-type: none"> 9. Basket Bracket- BC35051 10. Basket- S35003-PL 11. Tiller Swivel Assembly- S350109
Tiller Locking Knob- S31028 12. Steering Post- S350101 13. Tiller Bearing Set- S31021 14. Front Fork- S350105-HD 15. Steering Axle- S350102 16. Front Bumper- S350104-4HD 17. Connecting Rod (ea)- S350103 18. Seat Post- S31049 19. Seat Post Pin- S350169717 20. 17.5" Seat Assembly- S35045-HD
17.5" Back Cover- S35045BC-HD
17.5" Seat Cover- S35045SC-HD
20" Seat & Seat Plate- S35045-HD20 21. Armrest Complete Left- VENT-29
Armrest Complete Right- VENT-30 22. Arm Pad Left- VENT-31
Arm Pad Right- VENT-32 23. Seat Plate- S31039-HD 24. Front Frame- S350108 25. Rear Frame- S31130-4HD 26. Foot Shroud- S350107
Front Decoration Blue Left- S3501BL-1
Front Decoration Blue Right- S3501BR-1
Front Decoration Red Left- S3501RL-1
Front Decoration Red Right- S3501RR-1 27. Main Shroud- S350110
Main Decoration Blue Left- S3501BL
Main Decoration Blue Right- S3501BR
Main Decoration Red Left- S3501RL
Main Decoration Red Right- S3501RR | <ol style="list-style-type: none"> 28. Rear Shroud- S350111
Lifting Handle- S35011FH
Rear Decoration Blue Left- S35011B-L
Rear Decoration Blue Right- S35011B-R
Rear Decoration Red Left- S35011R-L
Rear Decoration Red Right- S35011R-R 29. Frame Contact Front- S31056HD
Frame Contact Rear- S31057HD 30. Frame Lock Handle- S350169716 31. Tension Knob- S350169715 32. Front Wheel Assembly- S350100-4HD
Front Wheel Bearing (pr)- S350100B-4HD 33. Rear Wheel Assembly- S31132B-HD 34. Controller (CTE)- S31008GT-HD
Controller (S-Drive)- C10-067-00200 *
* After Serial # (21)2W1801100007 35. Wire Harness (CTE)- S31048-HD
Wire Harness (S Drive)- S18-034-00400 36. Battery Pack Assembly- S350169707-HD
Battery Case Only- S350169706-HD
Battery (need 2)- LRM412123 37. Controller Cover- S350113 38. Anti-tip Wheel- S31026 39. Mud Flap Left- S35048MFL
Mud Flap Right- S35048MFR 40. Battery Charger- LRM302111 41. Motor-Transaxle Assembly- S31125HD
Motor-Brake- S31031-HD
Transaxle- S350112-HD
Transaxle Bracket- S350112BKT-HD
Free Wheel Knob- S31031K
Motor Axle Key- 29002W7011 42. Battery Contact- S31060 43. Bat. Box Charger Port- S350169706-HDCP 44. Armrest Knob- S31052 45. Front-Rear Axle Bolt- S350112AN |
|---|--|