

Titan P22 Front Wheel Drive Power Wheelchair Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Powerchair with Seating

Pan Seat

- TITAN18P22.....18" (W) x 18" (D)\$5,725.50
HCPCS: K0822, Includes interchangeable color shrouds in Blue & Red.
- TITAN20P22.....20" (W) x 18" (D)\$5,891.00
HCPCS: K0822, Includes interchangeable color shrouds in Blue & Red.

Standard Armrests

- ARM-L.....Included
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests
- ARM-R.....Included
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests

Controller Arm

- Adjustable Length.....Included
- AA4100 – Swing-Away Controller Arm (HCPCS E1028)\$105.00

Front Riggings

- Adjustable FootplateIncluded
- AE2500/LK3JELR – Elevating Legrest Kit\$207.48
HCPCS: E0990, Includes swing-away elevating legrest with attachment bracket

Batteries & Charger

- U1BATTERY (pair)Included
HCPCS: E2365, 12V x 35AH
- 3A Battery ChargerIncluded

Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978)Included
- SF8030 Crutch/Cane Holder (HCPCS E2207)\$120.12
- SF8010 Cylinder Tank Carrier (HCPCS E2208)\$156.52

Specifications

HCPCS

TITAN18P22	K0822
TITAN20P22	K0822

Performance

Maximum Weight Capacity	300 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	25"

Dimensions

Overall Length	40"
Overall Width	25"
Floor to Seat Height	20.5"–23"
Ground Clearance	2.5"

Power

Motor	Two 24V x 200W x 3900rpm
Controller	Dynamic Shark 40A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires

Casters	6" x 2" Solid
Drive Wheels	10" x 3" Flat Free

Weights

Base Weight	68 lbs.
Battery Weight	44 lbs.
Seat Weight	43 lbs.
Total Weight	155 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.