

Tamper-Proof Magnetic Pull Cord Alarm

HCPCS:



- When magnetic pull switch is removed, alarm is activated
- Pull switch can be activated from any direction
- Convenient On/Off switch for easy activation and deactivation
- Easily mounts on wheelchair or bed
- Cord adjusts from 28" 58"
- Alligator clip allows for attachment to patients clothes
- 97 103 dB
- One 9V battery included

PRODUCTS

Description

13603	Tamper-Proof Magnetic Pull Cord Alarm, EA
-------	---

ACCESSORIES

Description

136DC	ALARM AC ADAPTER, 1/EA, 25/CV, CS
13617	Cord,Magnet,Clip for 13603,1/ea, EA
13618	Replacement Plastic Carry,1/ea, EA