

## Lumbar Support General Use Wheelchair Back Cushion

HCPCS: E2611



FIGURE A



FIGURE B

- Memory-foam back cushion provides support, comfort and a tactile curve to improve posture
- Durable Masonite board prevents hammocking in sling-back wheelchair Easily removed as per patient preference
- Conforms to lumbar curve to support lower back
- Sealed liner protects foam core from moisture
- Nylon stretch cover is water resistant, non-slip vinyl back prevents the cushion from sliding and is waterproof
- Adjustable Straps easily attach to wheelchair
- Meets flammability standard TB117

### PRODUCTS

#### Description

14889	GeneralUseBackCushion,18"x17",1/cs
14906	GeneralUseBackCushion,16inx17in,1/cs
14920	General Use Back Cushion,20x17,1/ea
8033	Premier One FoamBackCushion 22x19 ,1/ea

### SPECIFICATIONS

#### Description

Bottom Fabric	Vinyl, waterproof
Top Fabric Material	Nylon stretch, water resistant
Warranty	18 Months Limited