

## Lumbar Support General Use Wheelchair Back Cushion

HCPCS: E2611



FIGURE A



FIGURE B

- Memory-foam back cushion provides support, comfort and a tactile curve to improve posture
- Durable Masonite board prevents hammocking in sling-back wheelchair Easily removed as per patient preference
- Conforms to lumbar curve to support lower back
- Sealed liner protects foam core from moisture
- Nylon stretch cover is water resistant, non-slip vinyl back prevents the cushion from sliding and is waterproof
- Adjustable Straps easily attach to wheelchair
- Meets flammability standard TB117

### PRODUCTS

### Description

|       |   |
|-------|---|
| 14889 | GeneralUseBackCushion,18"x17",1/cs      |
| 14906 | GeneralUseBackCushion,16inx17in,1/cs    |
| 14920 | General Use Back Cushion,20x17,1/ea     |
| 8033  | Premier One FoamBackCushion 22x19 ,1/ea |

### SPECIFICATIONS

### Description

|                     |                                |
|---------------------|--------------------------------|
| Bottom Fabric       | Vinyl, waterproof              |
| Top Fabric Material | Nylon stretch, water resistant |
| Warranty            | 18 Months Limited              |