

## Titanium Skin Protection and Positioning Gel/Foam Wheelchair Cushion

HCPCS: E2607




FIGURE A



FIGURE B

- The Titanium cushion is specifically designed to offer superior pressure redistribution and positioning
- High-resilient, high-density foam combined with a dual chamber gel bladder encourages proper patient posture through immersion
- Middle contour acts as a built-in leg abductor
- Low-shear, fluid-resistant stretch cover helps to maintain sensitive skin tissue
- Non-skid bottom offers additional stability
- Meets flammability standard TB117

### PRODUCTS

#### Description

PRODUCTS	Description
FPT-1	TitaniumViscoelasticfoam/GelCush16x16x3"
FPT-1818	Titanium Gel/Foam W/C cushion 18X18X3.5
FPT-2	TitaniumViscoelasticfoam/GelCush18x16x3"
FPT-2018	Titanium Gel/Foam W/C cushion 20X18X3.5
FPT-4	TitaniumViscoelasticfoam/GelCush22x18x3

### SPECIFICATIONS

#### Description

SPECIFICATIONS	Description
Cover Material	Fluid Resistant Stretch Top   Non-Skid Bottom
Warranty	1 Year Limited
Weight Capacity	Up to 20" 275 lbs.   22" and Up 500 lbs.