

## PreserveTech Universal Raised Toilet Seat



HCPCS: E0244



FIGURE A



FIGURE B



FIGURE C

- Accommodates three different seat heights (2", 4" and 6") by maneuvering just two components.
- Pivoting armrests can be raised out of the way or removed depending on your needs or bathroom environment.
- Fits comfortably on most standard and elongated toilets.
- Includes a lid with recessed grip for better hygiene and aesthetic and effortless opening.
- Contoured to provide a safe and comfortable elevation.
- Made from PP that is both lightweight and durable to support up to 300 lbs.
- Drive PreserveTech products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces.
- \*Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms.

### PRODUCTS

### Description

|              |                                       |
|--------------|---------------------------------------|
| RTL12C002-WH | PreserveTech Universal RTS, 1EA, 1/CV |
|--------------|---------------------------------------|

### SPECIFICATIONS

### Description

|                      |                  |
|----------------------|------------------|
| Seat Depth           | 19"              |
| Seat Width           | 16"              |
| Warranty             | Limited Lifetime |
| Weight Capacity      | 300 lbs          |
| Width Between Arms   | 18.5"            |
| Material Composition | Polypropylene    |