

AquaSense Bath Mat, Contoured with Temperature Indicator

HCPCS:



FIGURE A

- Built-in temperature indicator will show if bath water is too cold, too hot, or just right
- Accuracy of the indicator is within 1 degree of the specified zone
- Over 300 suction cups provide an active, secure grip
- Gentle, low profile textured surface for a secure grip where you stand
- Antimicrobial formulation helps control mildew and bacteria
- Accommodates specialized bath seating without compromising stability
- Contoured ends allow for better tub coverage without blocking the drain
- Gentle, low-profile surfaces provide a comfortable seating area

PRODUCTS

Description

785-530

AQUASENSE MAT W/TEMP.IND, 1EA,6CV,RET