

AquaSense Adjustable Bath Seat with Backrest



HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Both lightweight and durable the rust-resistant aluminum frame supports up to 300 lbs.
- Each leg includes a non-skid, non-slip tip offering security on wet surfaces.
- Simply and quickly adjust the height with the help of a push pin to achieve the ideal seated position. Seat height adjustments on this model range from 14" to 21".
- Drainage holes help to prevent water from pooling, making the bench more comfortable, less slippery and ensuring transfer from bathing to exiting the tub is more secure.
- Textured seat surface helps to prevent slipping.
- Push button, snap together design requires no tools and assembles in just minutes.
- Feel good knowing that plastics parts of this product are made using recycled materials.

PRODUCTS

	Description
770-510	AQS B.SEAT W/ BK WHITE, 1EA, 2CV, Ret.

SPECIFICATIONS

	Description
Waranty	Limited Lifetime
Weight Capacity	300 lbs