

AquaSense Adjustable Bath Seat without Backrest

HCPCS:



FIGURE A



FIGURE B

- Provides support for individuals who need to sit down while having a bath or a shower
- Textured non-slip surfaces with integrated drainage holes, as well as handles
- Corrosion resistant aluminum frame
- Adjustable legs raise or lower the seat height
- Rubber safety tips
- Tool-free installation
- Adjustable seat height from 14" to 21" (36 to 53 cm)
- AADL# B024

PRODUCTS

Description

770-500	AQS B.SEAT W/O BK WHITE, 1EA, 2CV, Ret.
---------	---

SPECIFICATIONS

Description

Weight Capacity	300 lbs
-----------------	---------