

Deluxe Aluminum Shower Bench without Back

HCPCS: E0245




FIGURE A



FIGURE B

- Tool-free assembly of seat and legs (Figure A)
- Drainage holes in bench reduce slipping
- Frame is lightweight, durable and corrosion proof
- New pinch-free cover allows for push pins to be depressed without pinching finger The audible "click" will be heard when secure (Figure B)
- Legs are height adjustable in 1" increments, and angled with suction style tips to provide additional stability
- Support collar prevents leg movement and reduces rattle

| PRODUCTS | Description |
|-----------|------------------------------|
| 12203KD-4 | Bath Seat w/o Back, KD, 4/cs |

| SPECIFICATIONS | Description |
|------------------------|----------------------------|
| Carton Shipping Weight | 7 lbs. 7 lbs. 21 lbs. |
| Outside Legs | 19.5" (W) x 19.5" (D) |
| Seat Dimensions | 19.25" (W) x 11.5" (D) |
| Seat Height | 15.5" - 19.5" |
| Warranty | Limited Lifetime |
| Weight | 5 lbs. 5 lbs. 4.5 lbs. |
| Weight Capacity | 300 lbs |

| ACCESSORIES | Description |
|-------------|---------------------------|
| 12002T | Tips for 12002,4/pkg., BX |