

## Transfer Tub Bench

HCPCS: E0247



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of back, legs and arm (Figure A - C)
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in  $\frac{1}{2}$ " increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

### PRODUCTS

#### Description

12011KD-1	Transfer Bench KD,Tool Free, 1/cs
12011KD-2	Transfer Bench KD Tool Free,2/cs

### SPECIFICATIONS

#### Description

Warranty	Limited Lifetime
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Weight Capacity	400 lbs.
Carton Shipping Weight	13 lbs.   13 lbs.   24 lbs.
HCPCS	E0247
Weight	10 lbs.

### ACCESSORIES

#### Description

12011SS	SuctionCups,Sm.12011,1pr/bx, BX
12011SL	SuctionCups,Lg.12011,1pr/bx, BX