

Bariatric Steel Forearm Crutch

HCPCS:



- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability
- Leg and forearm sections adjust independently for optimal sizing
- Vinyl hand grips are comfortable and durable
- Extra-large tips for added stability

PRODUCTS

Description

PRODUCTS	Description
10403HD	ForearmCrutch,Bariatric,Adult,1/cs

SPECIFICATIONS

Description

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight	5.7 lbs.
Weight Capacity	500 lbs