

## Mayo-Instrument Stand

HCPCS:



- Single post stand with "U" base design provides support and versatility
- Removable stainless steel tray measures 19" x 16 5"
- Tray height adjusts from 32" - 49" with lock
- Two, 3" casters provide additional mobility

PRODUCTS	Description
13035	Mayo-Inst.Stand,Single Post,1/cs
SPECIFICATIONS	Description
Base Dimensions	22" (W) x 18.5" (D)
Carton Shipping Weight	11 lbs.
Height	32" - 49"
Tray Dimensions	19" (W) x 13 " (D)
Warranty	Limited Lifetime
Weight	9 lbs.
ACCESSORIES	Description
13035T	Tray for 13035, 13045,13071,1/ea, EA