

Hydraulic, Deluxe Silver Vein Patient Lift

HCPCS: E0630



- Silver Vein steel construction
- 6 point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base

| PRODUCTS | Description |
|----------|---|
| 13023SV | Patient Lift,New Style,Silver Vein,1/cs |

| SPECIFICATIONS | Description |
|------------------------|---|
| Base Height | 6.5" (13023SV) |
| Base Width | 42" (open) / 22.5" (closed) |
| Boom Height | 17.5" - 72" (13023SV) |
| Carton Shipping Weight | 84 lbs. |
| Casters | 5" (13023SV) |
| Warranty | Limited Lifetime (frame), 1 Year Limited (pump) |
| Weight | 72 lbs. |
| Weight Capacity | 450 lbs |

| PARTS | Description |
|----------------|---|
| 13023SV01 | HangerAttachmentBolt13023SV,1ea, EA |
| 13023SV03 | MastComplete13023SV,1ea, EA |
| 13023SV04 | UpperMast13023SV,1ea, EA |
| 13023SV06 | PumpAttachmentBolt13023SV,1ea, EA |
| 13023SV07 | MastHingeBolt13023SV,1ea, EA |
| 13023SV08 | PushHandleGrip13023SV,1PR, PR |
| 13023SV09 | Pump13023SV,1ea, EA |
| 13023SV09-FPHH | Pump Handle Hardware, 13023SV, 1/ea, EA |
| 13023SV10 | SpreaderBarHandleGrip13023SV,1ea, EA |
| 13023SV11 | SpreaderBar13023SV,1ea, EA |
| 13023SV12 | ThumbScrew13023SV,1ea, EA |
| 13023SV13 | SpreaderBarReceiver13023SV,1ea, EA |
| 13023SV16 | RearCasterWithLock13023SV,1ea, EA |
| 13023SV18 | FrontCasterW/oLock13023SV,1ea, EA |
| 13023SV19 | BaseTipCover13023SV,1ea, EA |

PARTS**Description**

| | |
|-----------|-----------------------------|
| 13023SV20 | TubePlug13023SV,1ea, EA |
| 13023SV21 | CradlePivot13023SV,1ea, EA |
| 13023SV22 | PumpHandGrip13023SV,1ea, EA |
| 13023SV23 | MastTipCover13023SV,1ea, EA |
| 13052RL | ReleaseLever,13052,1/ea, EA |