

## LEVANTAR Floor Lift, Power Base 500 lbs, 227 kgs

HCPCS: E0635



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F



FIGURE G



FIGURE H

- Weight capacity of 500 lbs, 227 kgs
- Base option: Power
- Elegant, contemporary design with soft lines but still heavy duty
- Features you expect in a Premium lift
- Oversized sealed rear casters roll easily and lock securely
- Easily lifts resident from floor
- Redundant controls allow for use even if hand control is inoperable
- Smart & Safety Technology LCD display Visual battery fuel indicator on Hand Pendant and Lift Control Module Auditable low battery alert fuel indicator Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm Audible Alert if the lift detects weight overload Maintenance monitoring: by months or lift cycles Audible Alert that maintenance is due
- Visual battery fuel indicator on Hand Pendant and Lift Control Module
- Auditable low battery alert fuel indicator
- Lift safety stop if resident or caregiver accidentally gets in the way of the boom arm
- Audible Alert if the lift detects weight overload
- Maintenance monitoring: by months or lift cycles
- Audible Alert that maintenance is due
- NEW Six-point cradle with spring clips offer a more secure and wider foot print reducing sling squeeze, shallow height
- offers greater height range
- Can be charged on, or off the lift
- Includes 1 battery and external charger stand

### PRODUCTS

#### Description

FLP500	Levantar Floor lift, power base, 500 lbs
--------	--

### SPECIFICATIONS

#### Description

Battery Pack Weight	7.9 lbs (5.0A), 4.9 lbs (2.9A)
Battery Warranty	1 yr
Front Caster Diameter	3"
Gross Weight	156.6 lbs
Load Capacity Max.	500lbs
Net Weight	132.2 lbs
Overall Height	76.3"
Warranty on Battery	1 yr
Warranty on Electronics	5 yr

### ACCESSORIES

#### Description

FLSCALE	Patient Lift Scale 660lbs capacity, EA
---------	--