

## Bariatric Aluminum Folding Walker- Two Button

HCPCS:



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design accommodates individuals up to 600 lbs
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Vinyl-contoured hand grip

**PRODUCTS****Description**

10220-1	Walker,Oversized,Adt.,1/ea
10220-1WW	Walker, AdltHD, 2Btn w/Whls, 1/cs
10223-2	Walker,JuniorHvyDty.2Btn.,2/cs

**SPECIFICATIONS****Description**

Carton Shipping Weight	12   21   18 lbs.
Overall Width	25.25"   24.5"
Warranty	Limited Lifetime
Weight	9   8 lbs.
Weight Capacity	600 lbs
Width	20.5" inside hand grip   23" inside back legs