

Airgo One Button Folding Walker

HCPCS:



- Economical and lightweight walker
- Features plastic glide tips, 1" anodized aluminum tubing, and angled front legs
- The walker front bar also acts as a support rail for use over a commode bowl
- The one-button folding mechanism makes the walker easy to open and close
- AADL# C030
- ADP# MW1GM0010

PRODUCTS

Description

770-121	Airgo Walker, Wh+Ski, Adult, Silver, 1EA
770-126	Airgo Walker, Wh+Ski, Small, Silver, 1E
770-160	Airgo Walker, Tips, Adult, Silver, 1EA
770-165	Airgo F.Walker, Tips, Small, Silver, 1E

SPECIFICATIONS

Description

Weight Capacity	300 lbs
-----------------	---------