



## PreserveTech Deluxe Folding Two-Button Walker With 5" Wheels

HCPCS: E0143,E0143

- PreserveTech Deluxe Folding Walker combines strength and durability with the lightweight ease of use. The walker features a unique push button mechanism for easy folding. Its U-shaped frame offers greater clearance while its enhanced cross-frame design provides you with greater stability. The 5" wheels allow for navigation on many types of terrain, while the contoured hand grip and adjustable height make the walker comfortable to use
- Durable and lightweight design easily folds for storage
- Sturdy aluminum construction ensures maximum strength while remaining lightweight
- Easily fold and unfold the walker with the dual push-button design, which may be operated by fingers, palms or side of hand
- Unique U-shaped frame design allows you to bring the walker closer to assist in standing, as well as improved gait pattern when walking
- Soft, vinyl contoured handgrip
- Using the marked peg holes, easily adjust the height in 1" increments for a customized fit
- 5" wheels are ideal for maneuvering both indoors and outdoors
- Rear glide cap allows for walker to slide easily and smoothly over most surfaces
- Folds for easy transport and the sides are detachable for storing in tight spaces
- Quick and tool-free assembly
- Drive PreserveTech products are manufactured with an antimicrobial\* agent that helps to minimize the growth of stain and odor-causing bacteria and mold on treated surfaces
- \*Antimicrobial properties are built in to protect the products. These products do not protect users or others against bacteria, viruses, germs or other disease organisms

### PRODUCTS

|         | Description                   |
|---------|-------------------------------|
| 10210-1 | Walker Adlt 2but w/5 whl.1/ea |

### SPECIFICATIONS

|                 | Description      |
|-----------------|------------------|
| Base Dimensions | Depth: 18"       |
| Color           | Silver           |
| Handle Height   | 32" - 39"        |
| Material        | Aluminum         |
| Overall Width   | 25"              |
| Product Weight  | 7.5"             |
| Shipping Weight | 9.5 lbs          |
| Warranty        | Limited Lifetime |
| Weight Capacity | 350 lbs          |
| Wheel Size      | 5"               |

### ACCESSORIES

|        | Description                 |
|--------|-----------------------------|
| 10200B | Walker Basket,1/ea, RTL, EA |

## ACCESSORIES

## Description

|             |   |
|-------------|---|
| RTL10131    | Walker Caddy,1EA,3/CS, EA                   |
| 10124       | Walker Tray, Universal, 1/ea, EA            |
| 10125       | Walker Tray, 1/bx, RTL, BX                  |
| 10121       | Tennis Ball Glides,1pr/bx, RTL, BX          |
| 10110       | Glide Ski for Walker,1pr./bx, RTL, BX       |
| RTL10254BFL | Univ Mobility Bag, Med, Blue Flrl, 1EA, EA  |
| RTL10254GR  | Univ Mobility Bag, Med, Army Green, 1EA, EA |
| RTL10254TFL | Univ Mobility Bag, Med, Trop Flrl, 1EA, EA  |
| RTL10255NV  | Univ Mobility Bag, Small, Navy, 1EA, EA     |
| RTL10255PK  | Univ Mobility Bag, Small, Pink, 1EA, EA     |
| RTL10256NV  | Univ MobilityBag, Med Flat, Navy, 1EA, EA   |
| RTL10256PK  | Univ MobilityBag, Med Flat, Pink, 1EA, EA   |
| RTL10256TL  | Univ MobilityBag, Med Flat, Teal, 1EA, EA   |