

## Viper Plus Reclining Wheelchair

HCPCS: K0004, E1226-1, K0004 + E1226 + K0195 OR E0990



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Built-in seat rail extension and extendable upholstery easily adjusts seat depth from 16" or 18" (Patent pending) (Figure A - C)
- Attendant wheel locks in a convenient and easily accessible location
- Spreader bar acts as a back stabilizer and a convenient push handle for the attendant
- Offset rear wheel position adds additional stability and equally distributes the patients weight while in the various reclining positions
- Hydraulic reclining mechanism allows for adjustments up to 180° (Figure D - E)
- Durable, flame-retardant nylon upholstery resists mildew and bacteria
- Carbon steel frame with double-coated chrome for an attractive, chip-proof, maintainable finish
- Composite Mag-style wheels are lightweight and maintenance free Also features 8" front casters
- Nylon, padded armrests provide added patient comfort, and feature the flip-back style for easy transfers
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Standard dual axle provides easy transition of seat-to-floor height from adult to hemi-level
- Front caster forks are adjustable in three positions
- Comes standard with swing-away elevating leg rests
- Composite foot plates are durable and lightweight
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight

### PRODUCTS

### Description

PLA416RBDDA	ViperGT,16"Recl,FlpBackRemDsk,1/ea
PLA416RBDFA	ViperGT,16"Recl,FlpBackRemFull,1/ea
PLA418RBDDA	ViperGT,18"Recl,FlpBackRemDesk,1/ea
PLA418RBDFA	ViperGT,18"Recl,FlpBackRemFull,1/ea
PLA420RBDDA	ViperGT,20"Recl,FlpBackRemDesk,1/ea
PLA420RBDFA	ViperGT,20"Recl,FlpBackRemFull,1/ea

### SPECIFICATIONS

### Description

Back Height	23"   33"
Carton Shipping Weight	53 lbs   55 lbs   57 lbs
Height	50"
Length with Riggings	42"
Product Weight with 16 Seat	47 lbs
Product Weight with 20 Seat	51 lbs
Product Weight with 18 Seat	49 lbs
Seat Depth	16" - 18"
Seat To Floor Height	17.5"/19.5"

<b>SPECIFICATIONS</b>	<b>Description</b>
Warranty	Limited Lifetime
Weight Capacity	300 lbs
Width	24"   26.5"   28"
Width Closed	12.5"
Arm Rest Type (s)	Flip Back Detachable Desk Arm   Flip Back Detachable Full Arm
<b>ACCESSORIES</b>	<b>Description</b>
STDS6008-1	Carry Bag,W/C,Oxygen,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
STDS831	Heel Loops, 1pr, BX
WASL	Limb Support,Left,1/cs, CS
WASR	Limb Support, Right,1/cs, CS
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
STDS860	Forwrd Stablzrs/Univ.Fit ,1pr/bx, BX
<b>PARTS</b>	<b>Description</b>
STDS4S4826	Back Upholstery, 16", CS Recliner, 1/e, EA
STDS4S4804	HeadrestAssembly,20", CS Recliner, 1/ea, EA
STDS4S4801	Headrest Pillow, CS Recliner, 1/ea, EA
STDS4S4827	Back Upholstery, 18", CS Recliner, 1/e, EA
STDS4S4802	HeadrestAssembly, 16", CS Recliner, 1/, EA
STDS4S4808	Backrest Cane, CS Recliner, 1/ea, EA
STDS4S4828	Back Upholstery, 20", CS Recliner, 1/e, EA
TRHG	Hand Grips/Universal,1pr/bx, BX
STDS4Y4702L	ReleaseHndle w/Cable,SSP18RBDDA,LFT,1/ea, EA
STDS4Y4702R	ReleaseHndle w/Cable,SSP18RBDDA,RHT,1/ea, EA
STDS4Y4715	Anti Folding ADJ Knob,SSP18RBDDA,1/ea, EA
STDS4S4803	HeadrestAssembly, 18", CS Recliner, 1/, EA
STDS4S4809	SpreaderBar,16",CS,Sentra,Viper Rec,1/ea, EA
STDS4S4810	SpreaderBar,18",CS,Sentra,ViperRec,1/ea, EA
STDS3006	Cable,Back Release,Reclining w/c, 1/ea, EA
STDS4S4811	SpreaderBar,20",CS,Sentra,Viper,Rec,1/ea, EA
STDS4Y4704L	Piston,allCompnts,SSPRBDDA,Left,1/ea, EA
STDS3J2400	Handrim,Cruiser3,1/ea, EA
STDS004609-9	Hand Rim Screw Set of 9,1Bag, BAG
STDS4Y4704R	Piston,allCompnts,SSPRBDDA,Right,1/ea, EA