

### 3" Foam Retail Cushion

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- 3" foam cushion in attractive retail package with cutout to feel cover
- Comes with durable, machine washable, jacquard zippered cover
- Weight Capacity 250 lbs
- Made of soft, premium quality foam
- Versatile use for any location where seated comfort is needed

**PRODUCTS**

**Description**

RTL14910	Foam Cushion 18 x 16 x 3,1EA,12/CS
----------	------------------------------------

**SPECIFICATIONS**

**Description**

Cover Material	Jacquard Knit Cover
Warranty	1 Year Limited
Weight Capacity	250 lbs
Weight Capacity (13046   13046-HD)	250 lbs.