

Airgo One Button Folding Walker

HCPCS:



- Economical and lightweight walker
- Features plastic glide tips, 1" anodized aluminum tubing, and angled front legs
- The walker front bar also acts as a support rail for use over a commode bowl
- The one-button folding mechanism makes the walker easy to open and close
- AADL# C030
- ADP# MW1GM0010

| PRODUCTS | Description |
|-----------------|--|
| 770-121 | Airgo Walker, Wh+Ski, Adult, Silver, 1EA |
| 770-126 | Airgo Walker, Wh+Ski, Small, Silver, 1E |
| 770-165 | Airgo F.Walker, Tips, Small, Silver, 1E |
| SPECIFICATIONS | Description |
| Weight Capacity | 300 lbs |