

## Mobility Bags

HCPCS:



- Dual attachment options: back strap attachment with secure clip or adjustable hanging straps
- Securely attaches to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Four front exterior pockets and key ring zipper for added security of belongings
- Side mesh pocket holds water bottle
- Large, interior zip pocket
- Fully-lined pockets and interior
- Treated with water resistant technology
- Spot clean with wet cloth
- 
- 

### PRODUCTS

### Description

RTL10254BFL	Univ Mobility Bag, Med, Blue Flrl, 1EA
RTL10254GR	Univ Mobility Bag, Med, Army Green, 1EA
RTL10254TFL	Univ Mobility Bag, Med, Trop Flrl, 1EA

### SPECIFICATIONS

### Description

Retail Box Dimensions	13.8" L x 1.8" W x 13" H
Dimensions	11" H x 11" W x 3" D
Weight	.8 lbs.
Weight Capacity	9 lbs.