

## Mobility Bags

HCPCS:



- Dual attachment options: back strap attachment with secure clip or adjustable hanging straps
- Securely attaches to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Four front exterior pockets and key ring zipper for added security of belongings
- Side mesh pocket holds water bottle
- Large, interior zip pocket
- Fully-lined pockets and interior
- Treated with water resistant technology
- Spot clean with wet cloth
- 
- 

### PRODUCTS

### Description

|             |                                         |
|-------------|-----------------------------------------|
| RTL10254BFL | Univ Mobility Bag, Med, Blue Flrl, 1EA  |
| RTL10254GR  | Univ Mobility Bag, Med, Army Green, 1EA |
| RTL10254TFL | Univ Mobility Bag, Med, Trop Flrl, 1EA  |

### SPECIFICATIONS

### Description

|                       |                          |
|-----------------------|--------------------------|
| Retail Box Dimensions | 13.8" L x 1.8" W x 13" H |
| Dimensions            | 11" H x 11" W x 3" D     |
| Weight                | .8 lbs.                  |
| Weight Capacity       | 9 lbs.                   |