

Walker Accessory Bag

HCPCS:



Take your items with you wherever you go using the functional and versatile bag that easily secures to most walkers. The front and rear pockets are perfect for keys or a cell phone, and the large lined interior pocket with O-ring zipper closure holds all your essential items. The bag is constructed of lightweight nylon material that will remain sturdy even with extensive use.

- Adjustable bag loop fasteners attach to wide range of mobility aids
- Spacious interior pocket with zipper closure and key ring
- Three front pockets hold essentials such as keys, snacks or cell phone
- Two back pockets allow items to be accessed from users side
- Lined interior and pockets
- Treated with water resistant technology
- Spot clean with wet cloth

| PRODUCTS | Description |
|------------|---------------------------------------|
| RTL10256NV | Univ MobilityBag, Med Flat, Navy, 1EA |
| RTL10256PK | Univ MobilityBag, Med Flat, Pink, 1EA |

| SPECIFICATIONS | Description |
|-----------------------|-------------------------|
| Weight Capacity | 6 lbs. |
| Dimensions | 9" H x 16" W x 2.5" D |
| Retail Box Dimensions | 6.1" L x 1.4" W x 11" H |
| Weight | .5 lbs |