

Cast Protectors

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Keeps casts and bandages dry while bathing or showering
- Stretch tight fit eliminates the need for pumps, straps and adhesive tapes
- Reusable design folds compactly
- Latex Free

PRODUCTS

Description

RTLPC23401	Cast Protector, 20" Arm, 1 PC/RTL, 50/CV
RTLPC23402	Cast Protector, 24.5" Leg, 1 PC/RTL, 50/CV

SPECIFICATIONS

Description

Case Weight	6.6 lbs 9.15 lbs
Length	20" 24.5"
Weight	.11 lbs .16 lbs