

## Check and Go Oxygen Analyzer

HCPCS:



- 3-digit LCD display
- Ultra-light weight & stable
- Maintenance free
- Extended battery life
- Easy calibration
- Auto off

| PRODUCTS               | Description   |
|------------------------|---|
| 18580                  | O2 analyzer & sensor, disp 1 ea   |
| SPECIFICATIONS         | Description   |
| Accuracy and Linearity | 1% of full scale constant temperature, R.H and pressure when calibrated at full scale |
| Display Range          | 0.0 - 99.9% Oxygen  |
| Operating Temperature  | 60° F - 104° F (15° C - 40° C)  |
| Resolution             | 0.1% Oxygen   |
| Response Time          | 90% of final value < 15 seconds @ 77° F (25° C)                                       |
| Warm-up Time           | None  |
| Warranty               | Limited 2 Year  |
| Weight                 | Approx 3 oz   |