

Checker® Pediatric Compressor Nebulizer

HCPCS: E0570, A7003, A7005, A7015, A7013



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Child-friendly design to encourage therapy compliance
- Available in yellow or pink
- Easy-to-replace filter is located under the front hood
- Electrical cord located under rear hood
- Powerful piston pump is efficient, dependable and quiet
- Complete with mouthpiece, pediatric aerosol mask, tubing, two replacement filters and instructions
- Hood latch doubles as a cup holder
- Retail packaged

PRODUCTS

Description

18041-P	CheckerNebPinkReuseDispNebKit,1/ea,12/cv
---------	--

SPECIFICATIONS

Description

Electrical Requirements	115 VAC 60Hz
Liter Flow	8 lpm
Maximum Pressure	30 PSI
Method of Nebulization	Piston Compressor
Operating Pressure	16 PSI
Particle Size	0.5 µm to 5 µm
Unit Dimensions	11.8" (W) x 7.7" (H) x 4.3" (D)
Warranty	Limited 5 Year
Weight	3.7 lbs.

ACCESSORIES

Description

18090F	NebFilter 10/bag, BAG
18031	Nebulizer Carry Bag, 1/ea, EA
18254	Reusable Neb Kit 10/case, CS
MASK 001P	Pediatric Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit,Inc Mouth pc&Tub50/cs, CS