

## Deluxe Hip-High Chair



HCPCS:



FIGURE A

- Designed for post-hip surgery residents
- Allows for sitting without bending at the hip
- Armrests make sitting down or pushing up to get out from the chair safe and easy
- Attractive vinyl seat and backrest are soft, washable, comfortable and water resistant
- Constructed of solid maple hardwood not a plywood veneer
- The footrest has three height adjustments to provide maximum comfort and an ergonomic fit (Figure A)

PRODUCTS	Description
17100	Hip-High Chair, Cherry w/Bl, 1/c

  

SPECIFICATIONS	Description
Carton Shipping Weight	44 lbs.
Floor to Arm Height	31"
Floor to Footrest Height (adjustable)	5.5", 7.5", 9.5"
Floor to Seat Height	26"
Overall Dimensions	24" (W) x 40" (H) x 25" (D)
Warranty	1 Year on Upholstery, Limited 2 Year Warranty on Frame
Weight	33 lbs.
Weight Capacity	300 lbs
Width Between Arms	19"