

## Extra Soft Cozy Cannula

HCPCS: E0431



- Designed with curved prongs and an extremely soft cannula with a variety of tubing lengths to accommodate a wide number of patients' needs
- Compatible with supplemental oxygen sources delivering up to 15 LPM
- Latex free

PRODUCTS	Description
SOFT 200	Ad COZY cannula only CS/50
SOFT 204	COZYcannula w/4tub CS/50
SOFT 207	COZYcannula w/7tub CS/50
SOFT 207 P	Pedi COZYcannula w/7tub CS/50
SOFT 214	Ad Cozy cannula w/14' tub, 50/cs
SOFT 225	Ad COZYcannula w/25' tub 25/cs
SOFT 250	Ad COZY cannula w/50 tub CS/20